

Feature and Capabilities	Caristix Reader	Caristix Conformance
	Free	Paid
<b>Profile Library</b> <ul style="list-style-type: none"> <li>Browse trigger events, segments, fields in HL7 profiles created by reverse-engineering and editing.</li> </ul>	✓	✓
<b>Reverse-Engineering Profiles</b> <ul style="list-style-type: none"> <li>Load message logs to generate an HL7 profile.</li> <li>5 minutes gets you the HL7 structure and content of a sending system.</li> </ul>	✓	✓
<b>Profile Editing</b> <ul style="list-style-type: none"> <li>Edit profiles you build through reverse-engineering.</li> </ul>	✓	✓
<b>Segment Library</b> <ul style="list-style-type: none"> <li>View a list of segments, fields, and their attributes such as length and repeatability.</li> <li>Check for Z-segments and custom fields.</li> <li>One glance tells you if a system is populating a field.</li> </ul>	✓	✓
<b>Data Type Library</b> <ul style="list-style-type: none"> <li>View data type information.</li> <li>Check for custom data types and definitions that can trip up your integration.</li> </ul>	✓	✓
<b>Table Library</b> <ul style="list-style-type: none"> <li>View field content in tables.</li> <li>Access real-world values used by sending systems all in one place.</li> </ul>	✓	✓
<b>Layout Library</b> <ul style="list-style-type: none"> <li>Change your profile layout, add new columns.</li> <li>Build issue trackers. Add mapping information and business rules.</li> <li>One profile lets you track an entire interface build, without having to rely on Excel or Word for non-HL7 content.</li> </ul>		✓
<b>Gap Analysis</b> <ul style="list-style-type: none"> <li>Compare two profiles to surface the gaps you need to bridge in the interface.</li> <li>Structural and content gaps become your to-do list for configuring the interface.</li> </ul>		✓
<b>Synchronize Word output and Caristix profile</b> <ul style="list-style-type: none"> <li>Edit the Word document containing your profile, and sync changes back to the Profile Library.</li> <li>Perfect for analysts who prefer to work in Word, while preserving a canonical version for developers.</li> </ul>		✓

## About Caristix

The average US hospital runs up to 100 IT applications. Not a single one of them can share patient information out of the box. So hospitals and vendors turn to data interfaces – 50 to 100 of them in an average hospital. Each interface can take months of painstaking manual work to set up.

Caristix has developed a software suite to automate manual interface work. Our software reads HL7 data and outputs a list of interface requirements. As a result, Caristix software can reduce months of work to a few days. Reduce interface deployment time by 50%, reduce hospital testing time by 75%, and cut interface maintenance time by 90%.

## Download

Download a white paper on new HL7 integration capabilities enabled by Caristix technology at <http://promo.caristix.com/li-whitepaper-offer/>

## Follow Us

twitter.com/caristix  
linkedin.com/company/caristix  
caristix.com/blog

## Contact Us

Caristix  
1400 St-Jean-Baptiste Avenue, Suite 204  
Quebec City, QC, Canada G2E 5B7

1-877-872-0027

caristix.com  
info@caristix.com

Copyright © Caristix 2011. All rights reserved.